

Plastics as a science



Kartell
Labware division



Safety box Materiale: PC

Completamente autoclavabile, è utilizzato per il trasporto di provette o contenitori feci-urine ove si renda necessaria una tenuta stagna. È dotato di guarnizione a tenuta in silicone e 4 chiusure laterali. La confezione comprende un manico in acciaio inox per facilitare il trasporto e una scheda di utilizzo con istruzioni approvate dall'Organizzazione Mondiale di Sanità. Simbolo di rischio biologico serigrafato.

Safety box Material: PC

Completely autoclavable, it is used to carry tubes or urine-stool containers under safe conditions. With sylicon gasket and 4 side-locks. Provided with a stainless steel handle and user sheet with instructions approved by WHO. Bio-hazard symbol on the box.

Art.	Dimensioni / Dimensions
569	330x175x180 h
8569 Set 4 maniglie - Set of 4 sideloks	
8570 Guarnizioni in SI - Sylicon gasket	



Supporto per safety box Materiale: Acciaio inox

Progettato per 10 contenitori urine fino a 200 ml e 3 contenitori feci da 30 ml.

Safety box rack Material: Stainless steel

Designed for 10 urine containers up to 200 ml size and 3 stool containers 30 ml size.

Art.	
570	



Frigo Box Materiale: Polistirolo Antiurto

Contenitori frigo realizzati in materiale termoplastico antiurto con isolamento termico mediante camere isolanti in poliuretano; omologato per alimenti e quindi garantito completamente atossico. Sono disponibili delle Unità ghiaccio con 1 Kg di liquido, art. n° 3707. Adatto per contenere 2 safety box art. 569 per il trasporto esterno tra laboratori.

Frigo Box Material: High Impact Polystyrene

Low temperature containers made of high impact thermoplastic material, thermally isolated by means of a polyurethane isolating jacket; non-toxic and therefore 100% food-safe. A 1 Kg cooling block, art. 3707, is available. Suitable to contain 2 safety box art. 569 for outside transport between different laboratories.

Art.	Capacità lt. Capacity lt.	Dimensioni int. Dimensions int.	Dimensioni est. Dimensions ext.
3708	45	435 x 263 x 364h	523 x306 x403h

Kartell spa
Labware Division